



**Sizing line**  
for easy and accurate sizing

**Extra velcro**  
for adjustment to thin necks

**Closed/cell grey foam padding**  
is non/absorbent and stain-repellent

**Flip chin piece**  
preshaped to fit the chin  
May be flipped down  
to allow access for the  
patient's mouth

**Safety lock buttons**  
allow for maximum reliability

**Size marking**  
four standard adult sizes

**Large trach hole**  
is increased with sizing

**"Latches in"**  
provide for ratchet adjustment

**16 adjustable**  
adult settings

## Ambu<sup>®</sup> **Perfit ACE**<sup>™</sup> Adjustable Collar for Extrication

### **Perfection**

Ambu Perfit ACE is a one-piece design that adjusts to any adult/size neck with 16 different settings ranging from Neckless to Tall sizes.

### **Flexibility**

Adjustable to 16 precise settings, as well as four standard adult sizes.

The Ambu Perfit ACE incorporates latches that allow the rescuer to size the collar exactly to the victim's neck size. Two safety lock buttons hold the collar in place when adjusted. It folds completely flat and frees up valuable storage space. The collar is radiolucent, CT and MRI compatible.