

DATASHEET

The Central Venous Cannulation Simulator is designed for learning and practicing CVC techniques, and incorporates superficial anatomical details and functional accuracy to simulate the clinical experience.

Anatomical landmarks

The sternal notch, sternocleidomastoid, clavicle, and other supporting features are present. The sternal notch "sinks" 1 - 1.5 cm when palpated.

The neck of the simulator is positioned to the left, making landmark identity quite easy visually.

Internal features

Internal features include a replaceable muscle and bone section, subclavian, internal jugular, external jugular, and carotid.

Use of a Swan-Ganz catheter is possible and changing of the internal tubing is extremely easy.

Fluid supply stand not included.

Central Venous Cannulation



KEY BENEFITS

- Practice of central venous cannulation techniques
- Palpation exactly the same as on a real patient
- Accurate external and internal anatomical landmarks
- Neck is positioned to the left for easy landmark identification
- Use of a Swan-Ganz catheter possible

ORDERING SPECIFICATIONS

Item no.

Item no.	Central Venous Cannulation
LF01087U	incl. torso with muscle bone insert and skin, 2 fluid supply bags, 6 subclavian tubes, 6 internal jugular veins, 6 external jugular veins, 1 carotid pulse tubing, connector, 12 strapping ties, needle, blood, powder, hard carry case
LF01078U	Cannulation Skin/Vein Replacement

Ambu

Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com